

INSURANCE LIST

*BLUE CROSS BLUE SHIELD

*UNITED HEALTHCARE

*AARP

BRIDGEWAY

COVINGTON

MOUNTAIN STATES

SOUTHWEST SERVICES

USA PREFERRED

WORKMAN'S COMP PLANS

*BCBS PLANS EXCLUDING PREFIXES XBU, XBL, XBZ.

*UHC PLANS INCLUDE COMPASS PLANS & ALL SAVERS ALTERNATE FUNDING PLANS, EXCLUDES HUMANA AND PACIFICARE.

*AARP EXCLUDES GROUP: HCFAC9

INSURANCE LIST IS SUBJECT TO CHANGE WITHOUT NOTICE.

CONTACT US FOR INSURANCE VERIFICATION AT 520-290-5947.

FAX ORDERS TO OUR INSURANCE

INTAKE OFFICE: 520-437-0228


CREATING LASTING PARTNERSHIPS FOR THE CARE OF YOUR LOVED ONES AND EMPOWERING OUR COMMUNITY WITH QUALITY SERVICE AND PRODUCTS.



AMS RETAIL STORE
4638 E GRANT ROAD
TUCSON, AZ 85712
M-F: 9AM-5PM
SAT: 9AM-2PM



AMS MAIN OFFICE & WAREHOUSE
2100 N WILMOT, STE 111
TUCSON, AZ 85712
M-F: 9AM-5PM



AMS RESPIRATORY OFFICE
2100 N WILMOT, STE 110
TUCSON, AZ 85712
M-F: BY APPOINTMENT ONLY

520.296.5925 MAIN - 520.296.3620 FAX



ALLIED MEDCOR SERVICES, INC.
HOME MEDICAL EQUIPMENT &
OXYGEN THERAPY

SERVING THE TUCSON AND SURROUNDING
AREAS SINCE 1986.

PART OF THE LOCAL **FIRST** ARIZONA
ORGANIZATION.



PRESSURE RELIEF SOLUTIONS



FOR QUESTIONS OR TO PLACE AN ORDER CONTACT:

RANDY OWENS
MARKETING LIAISON
520.272.6002




PRICE \$142.98

MED-AIRE ALTERNATING PRESSURE PAD & PUMP SYSTEM
THE APP UTILIZES INFLATING AND DEFLATING AIR-FILLED BLADDERS TO AUTOMATICALLY CHANGE THE PRESSURE POINTS BENEATH THE PATIENT. THE 4 LPM PUMP PRODUCES CONSISTENT AIR FLOW AND PRESSURE OVER A 5 MINUTE CYCLE TIME THAT IS FIXED OR VARIABLE.

PREMIUM GUARD GEL OVERLAY
THE GUARD GEL OVERLAY ADDS 3.5" OF GEL FOAM TO THE MATTRESS. THE OVERLAY CONTAINS 15 GEL BLADDERS FOR THE PREVENTION AND TREATMENT OF PRESSURE ULCERS. VINYL BASE IS STAIN RESISTANT AND WATERPROOF. SMALLER VERSIONS ARE AVAILABLE FOR USE WITH WHEELCHAIRS.



PRICE \$204.00



PRICE \$580.00

ROHO DRY FLOATATION MATTRESS OVERLAY SYSTEM
COMPRISED OF SOFT, INTERCONNECTED AIR CELLS, THE DRY FLOATATION OVERLAY CONFORMS PRECISELY TO THE BODY'S CONTOURS, HELPING SOFT TISSUE TO MAINTAIN ITS SHAPE IN ORDER TO MAXIMIZE THE HEALING PROCESS AS WELL AS PREVENTION.

FOAM EGGCRATE MATTRESS PAD
CONVOLUTED SURFACE HELPS WITH WEIGHT DISTRIBUTION AND AIR CIRCULATION. IDEAL FOR PREVENTION AND TREATMENT OF DECUBITUS ULCERS. PROVIDES COMFORT AND SUPPORT FOR THOSE WHO ARE BED RIDDEN. MADE OF MEDICAL-GRADE FOAM FOR LONG-LASTING USE. SMALLER VERSIONS AVAILABLE FOR WHEELCHAIR USE.



PRICE \$45.99

WHAT SETS AMS APART FROM OTHER COMPANIES?

- HIGH DEDUCTIBLE, OR COPAYS, OR DON'T WANT TO GO THROUGH INSURANCE? WE HAVE THE LOWEST PRIVATE-PAY PRICES IN TOWN!
- NO LONG WAITS! - FREE SAME-DAY DELIVERY ON ORDERS \$50 OR MORE TO TUCSON METRO AREA
- OUR PROFESSIONAL TECHNICIANS SETUP & DEMONSTRATE OPERATION OF ALL OUR EQUIPMENT
 - KIND & COMPASSIONATE STAFF
- LOCAL-REPRESENTATION ONLY, NO CALL CENTERS
 - PAYMENT PLANS & OPTIONS AVAILABLE

CALL NOW FOR PRODUCT AVAILABILITY AND PRICING!



LIKE US ON

facebook



SCAN THIS QR CODE AND BROWSE OUR WEBSITE 24/7!



MED-AIRE 8" ALTERNATING PRESSURE & LOW AIR LOSS MATTRESS SYSTEM
ACTIVE MATTRESS REPLACEMENT SYSTEM PROVIDES BOTH ALTERNATING PRESSURE AND LOW AIR LOSS TO OPTIMIZE PRESSURE REDISTRIBUTION, SHEAR/FRICTION REDUCTION AND MICROCLIMATE CONTROL. DESIGNED TO PREVENT, TREAT AND HEAL PRESSURE ULCERS. CONTROL UNIT OFFERS 10 MINUTE CYCLES WHEN ALTERNATING AND CAN ALSO BE USED IN STATIC MODE. ADJUSTABLE PATIENT WEIGHT SETTINGS ALLOW FOR OPTIMAL IMMERSION AND PATIENT COMFORT & COMPLIANCE. SYSTEM INCLUDES 20 INDIVIDUAL AIR CELLS TO MAXIMIZE IMMERSION AND ENVELOPMENT.
PRICE \$804.00 / RENT: \$175.00/MO



PRICE \$59.99-99.99

GEL/FOAM WHEELCHAIR CUSHION
THE GEL/FOAM CUSHION HAS A HIGH-DENSITY FOAM SHELL SURROUNDING A UNIQUELY-DESIGNED GEL BLADDER. DUAL CHAMBERED GEL BLADDER PROVIDES PROPER IMMERSION AND ENVELOPMENT FOR PRESSURE REDISTRIBUTION. COVER IS WATER RESISTANT AND LOW SHEAR STRETCH TO REDUCE HEAT AND MOISTURE BUILD-UP.

ROHO WHEELCHAIR CUSHION
THE ROHO CUSHION FOR WHEELCHAIRS REDISTRIBUTES PRESSURE USING SEVERAL 3", INTERCONNECTED AIR CELLS THAT ARE EASY TO SET UP AND ADJUST. PROVIDES A COMFORTABLE, FUNCTIONAL AND STABLE SEATING ENVIRONMENT ALL WHILE HELPING TO PREVENT PRESSURE ULCERS.



PRICE \$204.00